



## One for Multi-use



### One for Multi-use

5L/min flow rate is suitable for supplying oxygen for anesthesia, anesthesia recovery, ICU, first aid etc.



### Cost-effective

Generate medical oxygen (USP93%) from the air. Cost effective and never run out of oxygen supply.



### Tri- filtration system

Ceramic fiber • Build-in filter • Humidifier bottle  
99.999% of bacteria, viruses and odor can be removed effectively.



### Well combined with other devices

Professional accessories to ensure that oxygen concentrator can be combined with devices like anesthesia machine, ICU house, mask and so on.



### High security


There are LOW O<sub>2</sub> PURITY ALARM and ABNORMAL ALARM to remind you to deal with the abnormal.




### Quiet

Quiet design is to offer you a silent and comfortable working environment.

ROC-5A Oxygen Concentrator is can be for scientific research and veterinary use. The high purity oxygen is separated from air by molecular sieve pressure swing adsorption (PSA) technology.



How the molecular sieve works



Oxygen concentration under the operational flow

## Specification

Item	Specification
Oxygen concentration	93%±3% (0-5LPM)
Power supply	230VAC~50Hz
Dimension	348*330*535mm
Oxygen tube	2m*2
Weight	18.5kg
Maximum outlet pressure	8.5psi
Average sound level	50.4dB(A)



# RWD

RWD Life Science Co.,Ltd

Add: 3F, Yuyang building, No.13 Hi Tech North Rd, Songping Mountain, Hi-tech Industrial Park North, Nanshan district, Shenzhen, Guangdong Province, China.  
Tel: (86)-755-86111286 ext 8283 Web: www.rwdstco.com E-mail: sales@rwdstco.com

RWD Life Science Inc.

Add: 6540 Lusk Blvd, Suite C161, San Diego, CA 92121  
Tel: (858)900-5879 E-mail: sales@rwdstco.com